

Eurys Series – Small Aperture (5mm)

Broadband Faraday Rotators and Isolators for Ti:Sapphire Lasers

EOT's *Eurys Series* Broadband Rotators rotate the plane of polarized light 90° at 800nm in the forward direction and 0° from 720–950nm in the reverse direction while maintaining the light's linear polarization. When placed between crossed polarizers, a broadband Faraday rotator becomes a broadband optical isolator. *Eurys Series* broadband optical isolators provide high transmission in the forward direction and strongly attenuates back-reflected light between 720–950nm in the reverse direction protecting Ti:Sapphire oscillators from the deleterious effects of back reflections and also eliminating preferential lasing at the lower gain wavelengths of Ti:Sapphire lasers. Utilizing optics with low refractive indices and short optical pathlengths minimizes pulse broadening due to dispersion in the optics associated with ultra-short laser pulses.



Benefits:

- Eliminate ASE from high-gain amplifiers that can cause parasitic or relaxation oscillations
- Prevent preferential lasing at low-gain wavelengths by providing broadband isolation

Features:

- 2° incident angle on polarizers for reduced re- turn losses
- Orthogonal isolated beams
- Precision pointing of isolated beams
- Precision mounting options
- Lockable port covers

Specifications^a:

	Rotator	Isolator ^b
Center Wavelength	800nm	800nm
Spectral Range	720–950nm	720–950nm
Polarizer Type	N/A	SF1 Glass
Transmission at 22°C	>98%	>92%
Isolation at 22°C	N/A	>33dB
Damage Threshold	>3.5J/cm ² at 10ns >1J/cm ² at 8ps	>3.5J/cm ² at 10ns >1J/cm ² at 8ps

^a Product specifications are subject to change.


^b Escape ports should be used if rejected light is >1W or 0.15J/cm² at 10ns or forward light is >25W. All stray beams should be properly terminated.

Note: All products are RoHS compliant.


Dispersion: Some pulse broadening does occur when using EOT's Broadband Isolators. The Sellmeiers Equation for TGG used in the broadband isolators is:

$$n^2 - 1 = \frac{E_d E_o}{E_o^2 - (hc / \lambda)^2}$$


where: E_o = 9.223eV and E_d = 25.208eV



Eurys Series Rotator Dimensions:




Eurys Series Isolator Dimensions:



Model Number:

A A	B B	C	D D D	E E E	F F F
Product Type	Aperture Size (mm)	Device type	Oper. Wavelength	Input Polarization	Output Polarization
BB	05	R-Rotator I-Isolator	800nm	000 045 090 135	000 045 090 135

Example: Description: 5mm isolator centered at 800nm; input horizontal, output vertical
 Model Number: BB-05-I-800-000-090



Input Polarization Reference

- Notes: 1) Light is rotated clockwise 90° from input to output for all catalog devices.
- 2) $\lambda/2$ waveplate available on output for arbitrary output polarizations, if required.

(07/25/2017)

Electro-Optics Technology, Inc.
 3340 Parkland Ct. Traverse City, MI 49686 USA
 (231)935-4044 | (800)697-6782 | Fax: (231)935-4046 | sales@eotech.com | www.eotech.com



Photonic Solutions Ltd Unit 2.2, Quantum Court, Research Avenue South,
 HWU Research Park, Edinburgh, EH14 4AP, UK, Tel: +44 (0)131 664 8122
 Email sales@photronicsolutions.co.uk Web www.photronicsolutions.co.uk